



Luminometer Tests

ATP Bioluminescence

Cloud Platform Automatic Syncing Input

CleanTrust™ ATP Surface Swab



Product No. YL0001




Applicable Instrument: YZ3000 Luminometer



Product No. YL0002

Applicable Instrument: YZ3000L & YZ3500 Luminometer

Suitable for SystemSURE PLUS / EnSURE Touch



Applicable Scope: Surface of the target subject

Sensitivity: 1×10^{-15} mol / 1×10^{-16} mol

Storage and Validity: 2~8°C, store away from light, valid for 12 months.

Low lot-to-lot variation, high accuracy and repeatability.

CleanTrust™ ATP Water Swab



WaterTotal

Product No. YL0003

Applicable Instrument: YZ3500 Luminometer

Suitable for SystemSURE PLUS / EnSURE Touch





WaterFree

Product No. YL0004

Applicable Instrument: YZ3500 Luminometer

Suitable for SystemSURE PLUS / EnSURE Touch



Applicable Scope: Water Samples

Sensitivity: 1×10^{-16} mol


Storage and Validity: 2~8°C, store away from light, valid for 12 months.

Total ATP = Free ATP + Microbial ATP

Microbial ATP = Total ATP - Free ATP

 Visual interpretation, no instrument needed.


CleanTrust™ Surface Protein Swab



Surface Protein Product No. YL0100

Sensitivity

Room temperature (20-25°C), 80µg after 1 minute,
50µg after 5 minutes, 20µg after 10 minutes.



Surface Protein Allergen Product No. YL0200

Sensitivity

5 µg after 15 minutes incubation at 37°C

Applicable Scope: Surface of the target subject

Storage and Validity: 2~8°C, store away from light, valid for 12 months.

☎ Tel: +86 532 5567 8026

✉ Email: info@yili-bio.com

📍 F4, Building D9, 1 Jinhui Road, National High-Tech
Development Zone, Qingdao, China.