



## Luminometer Tests

ATP Bioluminescence

Cloud Platform Automatic Syncing Input

### CleanTrust™ ATP Surface Swab



Product No. YL0001




Applicable Instrument: YZ3000 Luminometer



Product No. YL0002

Applicable Instrument: YZ3000L & YZ3500 Luminometer

Suitable for SystemSURE PLUS / EnSURE Touch



Applicable Scope: Surface of the target subject

Sensitivity:  $1 \times 10^{-15}$  mol /  $1 \times 10^{-16}$  mol

Storage and Validity: 2~8°C, store away from light, valid for 12 months.

Low lot-to-lot variation, high accuracy and repeatability.

### CleanTrust™ ATP Water Swab



**WaterTotal**

Product No. YL0003

Applicable Instrument: YZ3500 Luminometer

Suitable for SystemSURE PLUS / EnSURE Touch





**WaterFree**

Product No. YL0004

Applicable Instrument: YZ3500 Luminometer

Suitable for SystemSURE PLUS / EnSURE Touch



Applicable Scope: Water Samples

Sensitivity:  $1 \times 10^{-16}$  mol


Storage and Validity: 2~8°C, store away from light, valid for 12 months.

Total ATP = Free ATP + Microbial ATP

Microbial ATP = Total ATP - Free ATP

 Visual interpretation, no instrument needed.


## CleanTrust™ Surface Protein Swab



**Surface Protein** Product No. YL0100

**Sensitivity**

Room temperature (20-25°C), 80µg after 1 minute,  
50µg after 5 minutes, 20µg after 10 minutes.



**Surface Protein Allergen** Product No. YL0200

**Sensitivity**

5 µg after 15 minutes incubation at 37°C

Applicable Scope: Surface of the target subject

Storage and Validity: 2~8°C, store away from light, valid for 12 months.

Distributed in Australia by Essential Scientific

M: 0400 598 271

E: [info@essentialscientific.com.au](mailto:info@essentialscientific.com.au)

[www.essentialscientific.com.au](http://www.essentialscientific.com.au)

