


CleanTrust™ Sponge Stick

The sponge stick consists of a pre-moistened, antimicrobial-free sponge and a detachable plastic handle, all packaged in a sterile bag. Sample collection can be done without touching the sponge, and the handle can be quickly detached to prevent contamination. The product can be customized with a diluent, ideal for contactless sampling.




Benefits

- Long handle reduces contamination risk.
- For surface sampling and hard to reach areas.
- Detachable handle ensures hygiene, preventing direct contact.
- Optimized method enhances sampling efficiency.
- Ideal for large-area sampling.
- Simplified packaging cuts down on waste.
- Sponge maintains cell viability.
- Ensures stable, accurate sample analysis.



Procedure



Sponge Stick 40mm x 40mm Polyurethane Sponge

Product No.	Description	Specification
CT-DE2510S	D/E Neutralizing Broth, 10mL	20 pcs / bag
CT-BPW2510S	Buffered Peptone Water, 10mL	20 pcs / bag
CT-NB2510S	Neutralizing Buffer, 10mL	20 pcs / bag
CT-PBS2510S	Phosphate Buffered Saline, 10mL	20 pcs / bag

Tel: +86 532 5567 8026 Email: info@yili-bio.com

F4, Building D9, 1 Jinhui Road, National High-Tech Development Zone, Qingdao, China.